

Stirnleuchte Clar-F 55 Headlight Clar-F 55

Klassische, komfortable und hochwertige Stirnleuchte nach Clar mit Spiegel-Ø55mm.

Unsere Stirnleuchten zeichnen sich durch hohe Qualität, **hohe Lebensdauer** und ein optimales Preis/Leistungsverhältnis aus. Die **optimale Lichtausbeute** in Blickrichtung wird durch den speziell geformten Spiegel und eine **Silberbeschichtung** in Kombination mit der besonderen Wendel in der Glühlampe erreicht. Auch tiefe Körperhöhlen können so optimal ausgeleuchtet werden. Der **Lichtfleck** ist mit einem Stellrad **fokussierbar**.

Der **ergonomisch** geformte Stirnreifen ist gepolstert und kann stufenlos verstellt werden. Daher bleibt diese Stirnleuchte auch nach längerer Tragezeit **außerordentlich bequem**.

Für den Betrieb ist ein Transformator mit min. 1000 mA notwendig. Die Stirnleuchte ist mit folgenden Spannungen lieferbar:

4 V/500mA: Art.-Nr. 08-15104
6 V/1000mA: Art.-Nr. 08-15106
10 V/500mA: Art.-Nr. 08-15110

Ersatzteile:

Lampen: Art.-Nr. 99-01240 4V / 500mA
Art.-Nr. 99-433 6V / 1000mA
Art.-Nr. 99-01210 10V / 500mA

Spiegel: Art.-Nr. 08-154

Technische Angaben:

Gewicht: net. 220g / brut. 500g
Spiegel-Ø: ca. 55mm, Kabellänge: 3,2 m
Stecker: nach 1,6m Chinch-Stecker;
Bananenstecker Ø4mm am Ende

Stirnband: PE, weiß

Verpackung: Karton

Zubehör:

Tischtransformator 220V/0-10V 2A
Art.-Nr. 10-100

Steckertransformator 220V/6V, 1A
Art.-Nr. 10-10216

Schalter (zwischensteckbar)
Art.-Nr. 10-10002

Classical high quality Headlight by Clar with comfortable headband and mirror-Ø 55mm

Our Headlights are designed and manufactured for the highest standard of quality, **long lifetime** and competitive price. The light is **optimised** by a special lamp and focussed by the **silver coated mirror**,

which offers the best reflection ratio. Leading the light exactly to same direction than the eyes sight, the user is enabled to inspect all deep holes of the body. The mirror can be adjusted by two ball joints and the light spot is **adjustable** by an adjustment wheel. The **ergonomic** headband includes a bolster and can be adjusted stepless. Even after hours this

08-15106

headlight is still **highly comfortable** to wear. For operation a transformer with min. 1000mA is required (please order as accessory if required). The headlight can be delivered in following voltages:

4 V/500mA: Item code 08-15104
6 V/1000mA: Item code 08-15106
10 V/500mA: Item code 08-15110

Spare parts:

Lamps: Item code 99-01240 4V / 500mA
Item code 99-433 6V / 1000mA
Item code 99-01210 10V / 500mA
Mirror: Item code 08-154

Technical data:

Weight: net. 220g / brut. 500g
Mirror-Ø: 55mm, Length of cable: 3,2m
Plugs: Chinch plug after 1,6m,
Banana plugs Ø4mm at the end

Headband: PE, white

Packing: Carton

Accessories:

Desk transformer 220V/0-10V 2A
Item code 10-100

Wall plug transformer 220V/6V, 1A
Item code 10-10216

Switch Adapter
Item code 10-10002