

High-end-line

Extreme Powerful Light

► **08-815**

High performance, focusable headlight with High-Power-LED. Wireless, light weight, comfortable headlight with two interchangeable LiPo battery packs at the headband for working up to 7 hours and more. Adjustable light intensity by „Touch dimmer“



► **10-10815**

Charger for 2x LiPo battery



► **08-205-815**

LiPo battery pack, 2x included



Light weight and high performance

► **08-825**

High performance, focusable headlight with High-Power-LED; no batteries at the headband for lowest weight and highest comfort. LiPo battery pack with belt clip, adjustable light intensity by „Touch dimmer“



► **10-10780**

Charger for LiPo battery pack



► **08-205-825**

LiPo battery pack, 1x included

Ausgabe 11.2017 | Layout und Umsetzung: www.formfix.de

High power
LED-Stirnleuchten

LED headlights

www.faromed.de



Basic-line

Fantastic light

- ▶ **08-795-S**
LED head light, high precision light spot, wireless, basic model with **built in battery**, including charger 10-10795




- ▶ **08-795-W**
White headband

- ▶ **10-10795**
Charger



Easy and quickly battery change

- ▶ **08-797-S**
LED head light, high precision light spot, wireless, basic model with easily **interchangeable battery**, including charger 10-10795



Easy and fast battery change

- ▶ **08-205-815**
LiPo battery pack, 1 x included



- ▶ **10-10815**
Charger for 2x LiPo battery
Option; not included



Comfort-line

Touch dimmer, switch off automatic, robust headband for short time use (e.g. ENT)



Automatic On/Off Switch



*Touch dimmer
Easy hygienic cleaning, no knob*



Easy and fast battery change

- ▶ **08-796-S-1**
LED head light, best for short time use (e.g. ENT), high precision light spot, wireless, basic model with 1 x changeable battery, including charger 10-10795

Smooth Headband

- ▶ **08-806-S-1**
Same as 08-796-S-1 but smooth headband for long procedures, touch dimmer, switch off automatic



- ▶ **08-205-815**
LiPo battery pack, 1 x included



- ▶ **10-10815**
Charger for 2x LiPo battery
Option; not included

