

# Stirnleuchte Clar-F 55 Headlight Clar-F 55

**Klassische, komfortable und hochwertige Stirnleuchte nach Clar mit Spiegel-Ø55mm.**

Unsere Stirnleuchten zeichnen sich durch hohe Qualität, **hohe Lebensdauer** und ein optimales

Preis/Leistungsverhältnis aus. Die **optimale Lichtausbeute** in Blickrichtung wird durch den speziell geformten Spiegel und eine **Silberbeschichtung** in Kombination mit der besonderen Wendel in der Glühlampe erreicht. Auch tiefe Körperhöhlen können so optimal ausgeleuchtet werden. Der **Lichtfleck** ist mit einem Stellrad **fokussierbar**.

Der **ergonomisch** geformte Stirnreifen ist gepolstert und kann stufenlos verstellt werden. Daher bleibt diese Stirnleuchte auch nach

längerer Tragezeit **außerordentlich bequem**.

Für den Betrieb ist ein Transformator mit min. 1000 mA notwendig. Die Stirnleuchte ist mit folgenden Spannungen lieferbar:

4 V/500mA: Art.-Nr. 08-15104  
6 V/1000mA: Art.-Nr. 08-15106  
10 V/500mA: Art.-Nr. 08-15110

#### **Ersatzteile:**

Lampen: Art.-Nr. 99-01240 4V / 500mA  
Art.-Nr. 99-433 6V / 1000mA  
Art.-Nr. 99-01210 10V / 500mA

Spiegel: Art.-Nr. 08-154

#### **Technische Angaben:**

Gewicht: net. 220g / brut. 500g  
Spiegel-Ø: ca. 55mm, Kabellänge: 3,2 m  
Stecker: nach 1,6m Chinch-Stecker;  
Bananenstecker Ø4mm am Ende

Stirnband: PE, weiß

Verpackung: Karton

#### **Zubehör:**

**Tischtransformator** 220V/0-10V 2A

Art.-Nr. 10-100

**Steckertransformator** 220V/ 4V/6V, 1A

Art.-Nr. 10-10215

**Schalter (zwischensteckbar)**

Art.-Nr. 10-10002

08-15104  
08-15106  
08-15110



**Classical high quality Headlight by Clar with comfortable headband and mirror-Ø 55mm**

Our Headlights are designed and manufactured for the highest standard of quality, **long lifetime** and competitive price.

The light is **optimised** by a special lamp and focussed by the **silver coated mirror**, which offers the best reflection ratio. Leading the light exactly to same direction than the eyes sight, the user is enabled to inspect all deep holes of the body. The mirror can be adjusted by two ball joints and the light spot is **adjustable** by an adjustment wheel. The **ergonomic** headband includes a

bolster and can be adjusted stepless. Even after hours this headlight is still **highly comfortable** to wear. For operation a transformer with min. 1000mA is required (please order as accessory if required). The headlight can be delivered in following voltages:

4 V/500mA: Item code 08-15104  
6 V/1000mA: Item code 08-15106  
10 V/500mA: Item code 08-15110

#### **Spare parts:**

Lamps: Item code 99-01240 4V / 500mA  
Item code 99-433 6V / 1000mA  
Item code 99-01210 10V / 500mA  
Mirror: Item code 08-154

#### **Technical data:**

Weight: net. 220g / brut. 500g  
Mirror-Ø: 55mm, Length of cable: 3,2m  
Plugs: Chinch plug after 1,6m,  
Banana plugs Ø4mm at the end

Headband: PE, white

Packing: Carton

#### **Accessories:**

**Desk transformer** 220V/0-10V 2A

Item code 10-100

**Wall plug transformer** 220V/ 4V/6V, 1A

Item code 10-10215

**Switch Adapter**

Item code 10-10002