

Ersatzlampen passend für Nikon-Mikroskope
Replacement lamps fitting Nikon microscopes

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting Model	Spannung Voltage	Watt Wattage	Socket Base
99-72249	70249		6 V	20 W	P21D
99-72818	71818		10 V	2 A	P28D
99-72826	71826		24 V	150 W	G6.35
99-72168	72168		15 V	150 W	GZ 6.35
99-72223	76223		21 V	150 W	GX 5.3
99-72451	76451		20 V	150 W	SFcX12-4
99-72513	76513		12 V	100 W	G 6.35
99-72516	76516		12 V	50 W	G 6.35
99-72906	76906		6 V	15 W	Ba20d
99-72098	77098		6 V	0.04 A	Midget
99-72903	77903		6 V	15 W	P25D
99-72904	77904		6 V	30 W	BA15D
99-72910	77910		6 V	15 W	P25D
99-72153	78153		17 V	450 W	SK19/36
99-72564	78564		39...45 V	50 W	SFa6-2
99-72580	78580		-	200 W	SFc10-4
99-72590	78590		-	200 W	SFc10-4
99-72099	79099		6 V	20 W	G 4
99-72727	79727		12 V	5 W	Bi-Pin
99-72542	MXA-20-542		240 V	30 W	
99-72316	MXA-20-316		6 V	15 W	Ba15d
99-720SE		SE	220 V	20 W	Ba15d
99-Y5761		Eclipse E400	6 V	30 W	G 4
99-S126	500WS	Xenon Blitzlampe NF-505			

Ersatzlampen passend für American Optical Mikroskope
Replacement lamps fitting American Optical microscopes

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting Model	Spannung Voltage	Watt Wattage	Socket Base
99-78624	624		6 V	10 W	G4
99-78033	1033		6.5 V	2.75 A	P30d
99-78120	1120		6 V	20 W	G4
99-78121	1121		6 V	20 W	G4
99-78631	11581		6 V	4.5 A	P30d

Ersatzlampen passend für Bausch & Lomb Mikroskope
Replacement lamps fitting Bausch & Lomb microscopes

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting Model	Spannung Voltage	Watt Wattage	Socket Base
99-Y1051	71-71-55		2.7 V	0.4 A	special

Ersatzlampe für Fotometer / Analyzer
Replacement lamp for Fotometer / Analyzer

Artikel-Nr Code-No.	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Socket Base
99-S021	340845	Kodak / Johnsen & Johnsen	18 V	115 W	G 5.3
99-S131	2076896001	Roche	12 V	100 W	Litzen
99-Y1029	0013 031 029	Eppendorf		40 W	NOVAL EF80
99-Y1030	Na/10		15 V	15 W	Pico 9
99-Y1031		Sockeladapter von Pico 9 auf P28s			
99-Y1035	670-14602	Shimadzu SPD-10A			
99-S067	RF-530	Shimadzu Fluoreszenzdetektor			
99-Y1049	DM 2-30-5		6 V	15 W	Litzen
99-78840	641600/705-0840	Boehringer / Hitachi	12 V	20 W	Litzen
99-Y1036	F4500	Hitachi			

UV-Lampen für Fototherapie (UV und IR) und Entkeimung

UV-lamps for Foto Therapy (UV and IR) and Germicide

Artikel-Nr Code-No..	Vergleichstyp Reference	Für Modell Fitting model	Spannung Voltage	Watt Wattage	Sockel Base
99-UV001	F85/100W-PUVA	PUVA 100, 1000, 3001, 3002, 3003, 6001, 6002, 701 UV 1000K, 3002K, 3003K, 7001K, 8001K		100 W	G 13
99-UV002	F85/100W-UV6	UV 100, 1000, 3001, 3002, 6002		100 W	G 13
99-UV003	F85/100W-UV21	UV 1000K, 3002K, 3003K, 7001K, 8001K		100 W	G 13
99-UV004	F85/100W-01	UV 3003K, 7001K, 8001K		100 W	G 13
99-UV012	TL100W/01	Philips brand		100 W	R17d
99-UV005	F8 T5 PUVA	PUVA 200		8 W	G 13
99-UV007	UVA 40 W	PUVA 800, UV 800K		40 W	G 13
99-UV008	TL 20 W/12	UV 800, UV 800K		20 W	G 13
99-UV009	TL 100 W/10 R	PUVA 1000		100 W	G 13
99-UV010	FR90 T12	PUVA 4000, 4001		90 W	G13
99-UV011	FR74 T12	PUVA 6000		74 W	G13
99-UV049	F15 T8-PUVA	PUVA 180, PUVA 500		15 W	G 13
99-UV088	F15 T8 UV 21	UV 180		15 W	G 13
99-UV042	F20T12/D	white light		20 W	G 13
99-UV052	F20T12/D	blue light		20 W	G 13
99-UV047	F20T12/BB	super attiniche		20 W	G 13
99-UV064	G6T5			6 W	
99-UV046	G15T8 / TUV 15W LL			15 W	
99-UV045	G30T8			30 W	
99-UV101	HNS 15W OFR			15 W	
99-UV109	TL4W/08			4 W	
99-UV102	TL-D 18W/860	white light		18 W	
99-UV063	TL20W/12			20 W	G 13
99-UV019	TL20W/52	56018202 / 2 M 17415/17		20 W	G 13
99-UV015	TL20/33RS	2 M 17414		20 W	G 13
99-UV054	TL25W/09N			25 W	G 13
99-UV082	TL32W/840			32 W	G 13
99-UV072	TL-D 36W/54	TLD 36-1/54		36 W	G 13
99-UV084	TL-D 58W/865			58 W	G 13
99-UV110	Dulux S/E 5W/21/840				
99-UV066	Dulux S 9W/71	PL 9/52 / WY114		9 W	G 23
99-UV095	Dulux D 13W/21-840			13 W	
99-UV039	Dulux L 18W/71	2 M 21010		18 W	2G 11
99-UV041	Dulux L 18W/31	2 M 21009		18 W	2G 11
99-UV107	Dulux L 24W/840			24 W	
99-UV098	Q 255 Z4				
99-UV022	Q 402 Z7	56033065			R7s
99-UV027	Q 405 Z4	Psorilux 3060/3070			
99-UV059	Q 406 Z7	2 M 18785			R7s
99-UV090	Q 450 Z7				
99-UV028	Q 724 Z4	Psorilux 5050/5060			
99-UV056	UV 400F	Dr. Hönle		400 W	R7s
99-UV108	UV 400F/2	Dr. Hönle		400 W	R7s
99-S089	K14L	WEDECO 32986			
99-UV050	6600-0262-200	Ohmeda	20 V	140 W	GX5.3
99-UV051	6600-0680-200	Ohmeda	20 V	140 W	Litzen
99-UV036	Ultra-Vitalux 300W		230 V	300 W	E 27
99-UV061	Ultramed 250W		230 V	250 W	R7s
99-UV055	Ultramed 400W		230 V	420 W	R7s
99-UV037	Thera DL 150W		230 V	150 W	E 27
99-UV038	Philips IR R95 E		230 V		E 27
99-UV044	Philips IR 250R		230 V	250 W	E 27